

File formats for NGS data

sylvain.legrand@univ-lille.fr

Plan du cours

- **1** Sequence formats
- 2 Data quality
- 3 Alignment formats
- 4 Annotations formats
- 5 Graphical data visualization

1 – Sequence formats

Fastq format

4 lines per sequence

```
@SEQ_ID
GATTTGGGGTTCAAAGCAGTATCGATCAAATAGTAAATCCATTTGTTCAACTCACAGTTT
+
!''*((((***+))%%%++)(%%%%).1***-+*''))**55CCF>>>>>CCCCCCC65
```

- Line 1 : ID
- Line 2 : sequence
- Line 3: +, sometimes followed by the repetition of the sequence identifier
- Line 4: quality
- Format used for reads



Fastq fastq, ID line

@EAS139:136:FC706VJ:2:2104:15343:197393 1:Y:18:ATCACG

EAS139	the unique instrument name
136	the run id
FC706VJ	the flowcell id
2	flowcell lane
2104	tile number within the flowcell lane
15343	'x'-coordinate of the cluster within the tile
197393	'y'-coordinate of the cluster within the tile
1	the member of a pair, 1 or 2 (paired-end or mate-pair reads only)
Υ	Y if the read is filtered (did not pass), N otherwise
18	0 when none of the control bits are on, otherwise it is an even number
ATCACG	index sequence





Fastq format, quality

- Quality score= Phred Score
 - For each base, measure the probability that the assigned base is false
 - The score of a given base (Q) is given by the following equation: :

$$Q = -10log10(p)$$

p: Estimated probability that the given base is wrong So a high score indicates a lower probability of error.

```
Q10: 1 error per 10 bases \rightarrow 90% accuracy
```

Q20: 1 error per 100 bases \rightarrow 99% accuracy

Q30: 1 error per 1000 bases → 99.9% accuracy

• The score is encoded in ASCII to lighten the file: Illimina = Phred+33



Fasta format

 An ID line starting with « > » followed by the sequence with or without line breaks (example: every 80 bp)

```
>949344 pacid=16033748 polypeptide=917831 locus=949344 ID=949344.v1.107 annot-version=v1.0 ATGGCGACAAATAAATTTTCATCTTGTTTCTTATTTTCGTTAATGGTGTTTTTTTCCATCCTTTTGCCACTGATTTCAGG GCAAATGATACCATGTCTACTGGGGAAATGTAAAAACACAAGGACATGCAATGCAATGCATCTTGCAAATCTAGAGGATACAAAG GAGGGGCTTGTATAAGCATGGACGTTCGCTCAAAAACCGGTGCTTATTGCTGCAAAGTTAGATTTGAATAA >945418 pacid=16033749 polypeptide=913905 locus=945418 ID=945418.v1.107 annot-version=v1.0 ATGCCACCAAATATCTACAGACTCTCCTTCTTCTTATCCCTACTCTGTTTCTTCTTCATTCCCTGGTTCTTCTTCTCTCGA CGAACAAGGTCAAGCTCTCTTTTGGCATGGAAGTCTCAACTGAACATCTCCGGCGACGCTTTTTCCTCCTGGCACGTCGCCG ACACATCTCCCTGCAACTGGGTCGGCGTAAAATGTAACCGTAGAGGTGAAGTTTCGGAGATACAACTCAAAGGCATGGAC TTGCAAGGTTCTCTGCCGGTGACTAGTCTCCCGGAGCCTCAAGTCTCTTACTTCCCTCACTTTATCTTCACTCAATCTTAC CGGAGTAAATCCCCCAAGGAGATAGGAGACTTTATTGAGCTTGAATTACTCGATTTATCGGATAATTCTCTCAGGCGATA TCCCTTGGGAAATCTTCAGGCTCAAGAAACTCAAGACTCTGTTTTTTGAACACTAACAATCTCGAAGGTCGGATTCCGATG GAGATTCCGATGGGAAATCTTCAGGCTCAAGAAACTCAAGACTCTTTTTCGATAACAACTCAACAATCTCGAAGGTCGGATTCCGATG GAGATTTGGGAAATCTTCCGGAGGATTCCGATG
```

Used in particular for assembled sequences




2 – Data quality

QC software example: fastqc

http://www.bioinformatics.babraham.ac.uk/projects/fastqc/


Input format: fastq; output format: html report

Command line : fastqc file_in.fastq





Adapters trimming



TruSeq Adapter, Index 3

5' GATCGGAAGAGCACACGTCTGAACTCCAGTCACTTAGGCATCTCGTATGCCGTCTTCTGCTTG

TruSeq Adapter, Index 4

5' GATCGGAAGAGCACACGTCTGAACTCCAGTCACTGACCAATCTCGTATGCCGTCTTCTGCTTG

TruSeq Adapter, Index 5

5' GATCGGAAGAGCACACGTCTGAACTCCAGTCACACAGTGATCTCGTATGCCGTCTTCTGCTTG

TruSeq Adapter, Index 6

5' GATCGGAAGAGCACACGTCTGAACTCCAGTCACGCCAATATCTCGTATGCCGTCTTCTGCTTG

http://support.illumina.com/content/dam/illumina-

 $support/documents/documentation/chemistry_documentation/experiment-design/illumina-customer-sequence-letter.pdf$



Adapters trimming, software example: cutadapt

https://cutadapt.readthedocs.org/en/stable/

CP99_inf_1.fastq en cours... cutadapt version 1.0

Input format: fastq, Output format: fastq

Commande line example: cutadapt -b adaptateur1 -b adaptateur2 -0 10 file_in.fastq >file_out.fastq

Command line parameters: -e 0.1 -O 10 -a TGGAATTCTCGGGTGCCAAGG CP99_inf_1.fastq

```
Maximum error rate: 10.00%
 Processed reads: 12201246
  Trimmed reads: 564676 (4.6%)
 Too short reads: 0 ( 0.0% of processed reads)
  Too long reads: 0 ( 0.0% of processed reads)
    Total time: 489.68 s
  Time per read:
                   0.04 ms
=== Adapter 1 ===
Adapter 'TGGAATTCTCGGGTGCCAAGG', length 21, was trimmed 564676 times.
Histogram of adapter lengths
length count
       47388
10
11
       39176
12
       30270
13
       22786
14
       16764
15
       11535
16
       7717
17
       5232
18
       3710
19
       2830
20
       2352
21
       374916
```



Cleaning of the sequences according to their quality, software example: prinseq

http://prinseq.sourceforge.net/

Input file: fastq; output file: 2 fastq files, one with the poor quality sequences and another with the cleaned good quality sequences

```
Commande line: perl prinseq-lite.pl -fastq file_in.fastq -min_len 50 -min_qual_mean 25 -trim_qual_right 20 -ns_max_n 0 -noniupac
```

Cutadapt_CP99_miR_1.fastq en cours...

Input and filter stats:

Input sequences: 7,781,005 Input bases: 193,463,851 Input mean length: 24.86

Good sequences: 2,820,137 (36.24%)

Good bases: 55,524,434 Good mean length: 19.69

Bad sequences: 4,960,868 (63.76%)

Bad bases: 137,904,528 Bad mean length: 27.80

Sequences filtered by specified parameters:

trim_qual_right: 1746 min_len: 1026091 max_len: 3921798 min_qual_mean: 2526 ns max n: 8707



3 – Alignment formats

SAM/BAM formats

- https://samtools.github.io/hts-specs/SAMv1.pdf
- SAM: Sequence Alignment/Map format
- BAM: compressed version of the file
- TAB-delimited text format consisting of a header section, which is optional, and an alignment section
- Header lines start with '@', while alignment lines do not

```
@HD
        VN:1.0 SO:unsorted
@SQ
        SN:Chr1 LN:30427671
                                                             Header
@SQ
        SN:Chr2 LN:19698289
@SQ
        SN:Chr3 LN:23459830
@SQ
        SN:Chr4 LN:18585056
@SQ
        SN:Chr5 LN:26975502
@SQ
        SN:chloroplast
                        LN:154478
@SQ
        SN:mitochondria LN:366924
        "ID:BowtieVN:0.12.8CL:"bowtie -v 0 --all --best --strata -S Athaliana_TAIR10.fa ./SRR9602/Cutadapt_SRR9602_prinseq_good_e9wD.fastq
@PG
SRR960237.9
                                                                                           AGGCCTCTACGAATTCATGAT JIIJJJJHFHDHHFFFFFC@C
                                                                                                                                                     MD:Z:21 NM:i:0
                        Chr3
                                 12098843
                                                 255
                                                                                                                                            XA:i:0
SRR960237.20
                                 8658085 255
                                                 23M
                                                                                                                                                     MD:Z:23 NM:i:0
                0
                        Chr2
                                                                                   ACGGAATAATGTAAAACTGTACA
                                                                                                                   CCCFFFFHHHHHJJJJJJJJJJ
                                                                                                                                            XA:i:0
SRR960237.28
                16
                        chloroplast
                                         136311 255
                                                          20M
                                                                                           GAATTCACCGCCGTATGGCT JJJJJJGHHHHHFFFFFCCC
                                                                                                                                            XA:i:0
                                                                                                                                                     MD:Z:20 NM:i:0
SRR960237.28
                        chloroplast
                                         102319 255
                                                          20M
                                                                                           AGCCATACGGCGGTGAATTC CCCFFFFFHHHHHHGIJJJJJ
                                                                                                                                            XA:i:0
                                                                                                                                                     MD:Z:20 NM:i:0
                                                             Alignment
```



Header in SAM/BAM

 Each header line begins with the character '@' followed by one of the two-letter header record type codes

Examples:

- @HD: VN: version of the file; SO: sorting order of alignments
- @SQ: Reference sequence dictionary. The order of @SQ lines defines the alignment sorting order
- @PG (Program): ID; VN: version; CL: command line
- Others fields possible...



Alignment in SAM/BAM

Each alignment line has 11 mandatory fields

Col	Field	Type	Regexp/Range	Brief description
1	QNAME	String	[!-?A-~]{1,254}	Query template NAME
2	FLAG	Int	$[0, 2^{16} - 1]$	bitwise FLAG
3	RNAME	String	$*[:rname:^*=][:rname:]*$	Reference sequence NAME ¹¹
4	POS	Int	$[0,2^{31}-1]$	1-based leftmost mapping POSition
5	MAPQ	Int	$[0,2^8-1]$	MAPping Quality
6	CIGAR	String	* ([0-9]+[MIDNSHPX=])+	CIGAR string
7	RNEXT	String	* = [:rname:^*=][:rname:]*	Reference name of the mate/next read
8	PNEXT	Int	$[0,2^{31}-1]$	Position of the mate/next read
9	TLEN	Int	$[-2^{31}+1, 2^{31}-1]$	observed Template LENgth
10	SEQ	String	* [A-Za-z=.]+	segment SEQuence
11	QUAL	String	[!-~]+	ASCII of Phred-scaled base QUALity+33



SAM/BAM handling

- samtools Utilities for the Sequence Alignment/Map (SAM) format
- http://www.htslib.org/doc/samtools.html
- Samtools is a set of utilities that manipulate alignments in the SAM and BAM formats.
- It converts between the formats, does sorting, merging and indexing, and can retrieve reads in any regions swiftly.



4 – Annotations format

BED format

- Browser Extensible Data
- Utilisé pour stocker des régions génomiques sous forme de coordonnées ainsi que les annotations associées
- The first three required BED fields are:
 - chrom: name of the chromosome (e.g. chr3, chrY, chr2_random) or scaffold (e.g. scaffold10671).
 - chromStart: starting position of the feature in the chromosome or scaffold. The first base in a chromosome is numbered 0.
 - chromEnd: ending position of the feature in the chromosome or scaffold (exclusive)
- Note: The chromEnd base is not included in the display of the feature, however, the number in position format will be represented. For example, the first 100 bases of chromosome 1 are defined as chrom=1, chromStart=0, chromEnd=100, and span the bases numbered 0-99 in our software (not 0-100), but will represent the position notation chr1:1-100



BED format

• 9 additional optional BED fields can be used (name, score, strand...)

chr7	127471196	127472363	Pos1	0	+
chr7	127472363	127473530	Pos2	0	+
chr7	127473530	127474697	Pos3	0	+
chr7	127474697	127475864	Pos4	0	+
chr7	127475864	127477031	Neg1	0	-
chr7	127477031	127478198	Neg2	0	-
chr7	127478198	127479365	Neg3	0	-
chr7	127479365	127480532	Pos5	0	+
chr7	127480532	127481699	Neg4	0	-



GFF format

- GFF (General Feature Format)
- GFF lines have nine required fields that must be tab-separated
 - seqname: The name of the sequence (chromosome or scaffold...)
 - source: The program that generated this feature.
 - feature: The name of this type of feature (ex: "CDS", "start_codon", ...)
 - start: The starting position of the feature in the sequence. The first base is numbered 1
 - end: The ending position of the feature (inclusive).
 - score: A score between 0 and 1000
 - strand: Valid entries include "+", "-", or "."
 - frame: If the feature is a coding exon, frame should be a number between 0-2 that represents the reading frame of the first base.
 - group: All lines with the same group are linked together into a single item.
- If a field is empty, enter "."



GFF format

##gff-version 3								
##annot-version v1.0								
##species Arabidopsis lyrata								
scaffold_1	phytozomev11	gene	29681	31614		+		ID=311229.v1;Name=311229
scaffold_1	phytozomev11	mRNA	29681	31614		+		ID=311229.v1.107;Name=311229;pacid=16057706;longest=1;Parent=311229.v1
scaffold_1	phytozomev11	exon	29681	29746		+		ID=311229.v1.107.exon.1;Parent=311229.v1.107;pacid=16057706
scaffold_1	phytozomev11	CDS	29681	29746		+	0	ID=311229.v1.107.CDS.1;Parent=311229.v1.107;pacid=16057706
scaffold_1	phytozomev11	exon	29864	29938		+		ID=311229.v1.107.exon.2;Parent=311229.v1.107;pacid=16057706
scaffold_1	phytozomev11	CDS	29864	29938		+	0	ID=311229.v1.107.CDS.2;Parent=311229.v1.107;pacid=16057706
scaffold_1	phytozomev11	exon	30016	30117		+		ID=311229.v1.107.exon.3;Parent=311229.v1.107;pacid=16057706
scaffold_1	phytozomev11	CDS	30016	30117		+	0	ID=311229.v1.107.CDS.3;Parent=311229.v1.107;pacid=16057706
scaffold_1	phytozomev11	exon	30232	30274		+		ID=311229.v1.107.exon.4;Parent=311229.v1.107;pacid=16057706
scaffold_1	phytozomev11	CDS	30232	30274		+	0	ID=311229.v1.107.CDS.4;Parent=311229.v1.107;pacid=16057706
scaffold_1	phytozomev11	exon	30444	30448		+		ID=311229.v1.107.exon.5;Parent=311229.v1.107;pacid=16057706
scaffold_1	phytozomev11	CDS	30444	30448		+	0	ID=311229.v1.107.CDS.5;Parent=311229.v1.107;pacid=16057706
scaffold_1	phytozomev11	exon	30605	30685		+		ID=311229.v1.107.exon.6;Parent=311229.v1.107;pacid=16057706
assets 1 d 1		CDC	COCOE	2000			0	TD 2112201 107 CDC C. Dannet 2112201 107. page 4 10057700




Handling of BED and GFF files

- Bedtools: A swiss-army knife of tools that allow a fast and flexible way of comparing large datasets of genomic features
- https://bedtools.readthedocs.io/en/latest/index.html
- Bedtools allows one to intersect, merge, count, complement, and shuffle genomic intervals from multiple files in widely-used genomic file formats such as BAM, BED, GFF/GTF, VCF
- The full list of bedtools sub-commands: annotate, bamtobed, bamtofastq, bed12tobed6, bedpetobam, bedtobam, closest, cluster, complement, coverage, expand, flank, fisher, genomecov, getfasta, groupby, igv, intersect, jaccard, links, makewindows, map, maskfasta, merge, multicov, multiinter, nuc, overlap, pairtobed, pairtopair, random, reldist, shift, shuffle, slop, sort, subtract, tag, unionbedg, window



Handling of BED and GFF files

Examples:






5 - Graphical data visualization - example: IGV

- The Integrative Genomics Viewer (IGV) is a high-performance, easy-touse, interactive tool for the visual exploration of genomic data
- https://software.broadinstitute.org/software/igv/
- Formats
 - BAM BED
 - BEDPE
 - BedGraph
 - bigBed
 - bigWig
 - Birdsuite Files
 - broadPeak
 - CBS
 - Chemical Reactivity Probing Profiles
 - chrom.sizes
 - CN
 - Custom File Formats
 - Cytoband
 - FASTA
 - GCT
 - CRAM
 - GFF/GTF

- Goby
- GWAS
- IGV
- LOH
- MAF (Multiple Alignment Format)
- MAF (Mutation Annotation Format)
- Merged BAM File
- MUT
- narrowPeak
- PSL
- RES
- RNA Secondary Structure Formats
- SAM
- Sample Info (Attributes) file
- SEG
- TDF
- Track Line
- Type Line
- VCF
- WIG



IGV







Sylvain Legrand
Maître de Conférences
UMR CNRS 8198 EVO-ECO-PALEO
Evolution, Ecologie et Paléontologie
Université de Lille - Faculté des Sciences et Technologies
Bât SN2, bureau 208 - 59655 Villeneuve d'Ascq

sylvain.legrand@univ-lille.fr | http://eep.univ-lille.fr/ Tél. +33 (0)3 20 43 40 16